User instructions

Carrara Volumia

0344

Dental ceramic system for over-pressing metal sub-structures for crowns and bridges.

Rx only
ISO 9693 + ISO 6872

Elephant Dental B.V.
Verlengde Lageweg 10
1628 PM Hoorn, The Netherlands
Tel: +31 229 25 90 00
Fax: +31 229 25 90 99
Email: info@elephant.nl
www.elephant-dental.com

Made in Holland

Elephant Hoom Holland

1. General information

1.1 Product description

Carrara Volumia is a new system component in the Carrara Original System that offers an interesting and economical alternative to the dental ceramic layering technique. The Carrara Volumia technique is designed to be the "full contour technique". The combination of the clear-colored *Linear* and the press pellets which transparency harmonizes between the dentin and the enamel enables the technician to make a reliable and simple reconstruction of the 16 classical "V" classes. Carrara Volumia can also be used in the Colback-technique.

1.2 Indication

Dental ceramic for over-pressing metal sub-structures for crowns and bridges made is:

- Carrara P&F
- Care E
- Vi-Comp LFC

Maximum wax weight (weight without sub-structure) of 1.4 g with a minimum strength of 0.8 mm (wax-ep).
Carrara P&F: The restoration should never exceed 4 units, with a maximum of one pontic.

1.3 Contra-indication

The over-pressing in combination with:

- Zirconium oxide sub-structures
- Aluminium oxide sub-structures
- Glass infiltrated aluminium oxide substructures
- All other alloys not listed under 'Indication'

Bruxism and other parafunctions

If patients are known to be allergic to one of the components, the material should not be applied.

1.4 Reuse

Carrara Volumia cannot be reused.

1.5 Precautions

Consult the Material Safety Data Sheet (MSDS) for more information.

2. Preparation

2.1 Preparation
To ensure a successful restoration using the Carrara Volumia system, please follow the general guidelines for preparation and layering thicknesses, as also described for full-ceramic restorations.

Preparation form	Circular preparation	Angle	Incisal/occlusal
Shoulder	1.5 mm (+/- 0.3 mm)	110°-130°	1.5 mm-2.0 mm
Chamfer	1mm (+/- 0.5 mm)	110°-130°	1.5 mm-2.0 mm
Knife-edge	-	110°-130°	1.5 mm-2.0 mm

¹ In occasion of pressing a ceramic glaze, the metal edge can be reduced with maximal 1mm.

3. Processing

3.1 Preparation of metal framework

- The anatomical shape of the restoration should be considered when designing the sub-structure, as is the case with the standard dental ceramic layering technique.
- To stabilize the metal sub-structure in the muffle after burnout, you must apply a retention bar (T-shaped lingual, palatal).
- Instead of oxidizing the alloy, a liner wash-bake must be used for all three alloys.
- Apply the Carrara Volumia Liner in a thin coat and fire according the firing schedule below (see 3.2).
- Afterward apply a modeling layer of Carrara Volumia Liner (see the user instructions for the full contour technique and the layering technique in the Carrara Interaction manual).
- If required, a second coat can be applied again afterwards.

It is also possible to use hollow pellets. It is recommended to fill the hollow pellets with Liner first, before the final ceramic layers are applied.

Preheat	Start temp.	Rate of heat increase	Final temp.	Vac.	Hold time
4 min.	400°C	60°C/min.	915°C	1 min.	2 min.*

¹ 1 minute with vacuum, 1 minute without vacuum.

3.3 Liner (paste opaque) firing chart

Start temp.	Rate of heat increase	Final temp.	Vac.	Hold time	
8 min.	400°C	60°C/min.	915°C	1 min.	2 min.*

¹ 1 minute with vacuum, 1 minute without vacuum.

4.1 Wax-ep

- Before starting the wax-ep, the restoration has to be weighed (incl. Liner) to determine the exact weight of the wax.
- Prepare the required wax model 'full contour' by using a fully residue-free combustible modelling wax suitable for press-ceramic systems.
- Keep in mind that the minimum thickness of the wax model should be no less than 0.8 mm.
- With the Carrara Volumia technique, slope and occlusion of the wax-ep can be perfectly modeled and checked in an articulator.

4.2 Layer thickness

Layer thickness sub-structure: min. 0.4 mm after finishing
To ensure the correct aesthetic phase keep to equally divided layer thickness among the bridge units.

4.3 Spraying

- It is important to keep the following guidelines in mind when spraying the press-ceramic:
 - When spraying with only one wax-ep procedure, the diameter of the wax wire should be at least 2.5 mm Ø using a fully residue-free combustible wax suitable for press-ceramic systems.
 - Do not exceed the length of max. 6 mm.
 - Change nozzle, attach press channel at every 30° arc (Ø 3 mm) to distribute the pressure evenly during pressing.
 - Preserve the occlusal contacts (I,II buccal, I,II lingual).
 - The muffle furnace must be kept completely clean and free of wax or other residues.
 - Check the correct position of the object (see picture 1&2).
 - Avoid sharp edges (see picture 1 & 2).
 - Ensure that all fractures blend (see picture illustration).
- Determine the final wax weight, including the pressure channel(s).
- The spraying angle from the wax object to the muffle furnace must be 45° (see picture 1&2).

Start temp.	Vacuum Level	Rate of heat increase	Press temp.	Hold time	Press time	Pressure
700°C	50 HPa	60°C/min.	940°C	20 min.	10 min.	2.7 bar

4.4 Wax weight chart including press channel	
1 pellet	2 pellets
Wax weight 150µg	up to 0.7 g
Number of objects	0.7 g – 1.4 g
Do not exceed the max. wax-weight.	Single crowns Bridges with max. 3 pontic or multiple crowns*

5.1 Investing

For investing, use Carrara Universal Denture Investment (Art. no. 232978050) and follow the user instructions. (Recommended liquid concentration: 50% Carrara Universal DL Liquid – 50% distilled water)

5.2 Burnout program (preheating)

Start temp.	Rate of heat increase	Final temp.	Hold time
600°C	10°C/min.	900°C	30 min.

- Please follow the manufacturer's instructions for your specific furnace.
- Minimal burnout time for the Carrara Volumia system: 30 minutes.
- Do not exceed the maximum burnout time of 2 hours!

- Note:
- For each additional muffle, the minimal hold-time should be increased by 10 minutes.
- The maximum hold-time should not exceed 2 hours!

6.1 Pressing

- Take the muffle out of the furnace immediately after preheating.
- Insert the required quantity of press pellets (see chart 4.4) and the Carrara Press Plunger (Art. no. 004100580) with the smooth side facing down and place the muffle in the furnace.

Colour-schedule for Full-Contour only

Indication	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Liner	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Transpa Dentite Pellet	T01	T01	T02	T02	T03	T04	T04	T05	T05	T02	T02	T03	T03	T01	T04	T04
Carrara Shade point	1	2	2	2	2	4	5	4	4	5	6	7	7	3	8	9

- Note:
- Please be aware that after inserting the first press pellet, it is advised to start by positioning the Carrara Press Plunger to ensure that the press pellet is correctly inserted.
- Then insert the second press pellet and secure it with the Carrara Press Plunger, as before.

6.2 Pressing parameters StrataPress

Start temp.	Rate of heat increase	Final temp.	Hold Time	Press time	Vac. Start	Vac. level	Vac. level	Pressure
700°C	60°C/min.	940°C	20 min.	10 min.	700°C	101%*	940°C	5 bar

* Continuous maximum vacuum

After pressing:

- Take the muffle out of the furnace and let it cool to room temperature.
- Do not force the cooling process (at pressure)
- After cool-down, direct immediately.

Note:

Press parameters for other furnace manufacturers see 11.

7.1 Divesting

- Place the height of the cylinder on the outside of the muffle with the help of another press plunger.
- Make an incision at the marked height using a separating disk.
- Now you can safely separate the ceramic restoration from the plunger.

7.2 Sandblast

- Direct the object by sandblasting with 50 µm glass beads (i.e. Elephant Bright Blast 50, Art. no. 255052310) at a max. pressure of 1.5 bar.

8.1 Finishing

- Carefully cut off the press channels with a diamond disc.
- To avoid local overheating of the ceramic, make sure that neither the working pressure nor the rotational speed (max. 8.000 rpm) of the disc is too high.
- The final trimming and shaping of the ceramic restoration can be done using sharp and fine-grained diamond cutters (max. 15.000 rpm).
- The entire surface of the restoration must be reworked.

8.2 Glaze firing

- Clean the restoration thoroughly with a steam cleaner before you start with the glaze firing.
- Use the Carrara Shade points (Art. no. 435003401 – ...409) to improve the natural appearance of the restoration and to copy other aesthetics. With these 9 shades, it is also possible to optimize the diameter of the 16 V classes. This can also be done using the Carrara Paint system (Art. no. 435095070) or the Classic Stains (Art. no. 49529510).

8.3 Shade point & Glaze firing program for StrataPlus / StrataPress

Dry	Start temp.	Vacuum Start	Rate of heat increase	Final temp.	Vacuum Off	Hold time
4 min.	450°C	450°C	60°C/min.	830°C	830°C	1-2 min. (without Vac.)

Glaze firing programs from other furnace manufacturers can be found in the extensive Carrara Interaction manual.

Note:

- In combination with Care E, cool down each firing cycle for 5 minutes.
- Carrara P&F & Vi-Comp LFC do not require any cooling.

9. "Carback-techniek"

- If the anatomical shape of the restoration (and/or other affect-materials) is desired, the required space can be provided if the dentin-cure is shaped and pressed in an anatomical shape instead of pressing a "full contour" crown.
- In the case of the part anatomized wax up (denture shade), make sure that the wax-ep has a minimal thickness of 0.8 mm and the wax weight does not exceed 1.4 g.
- All further procedures correspond to those in the user instructions for the full contour technique and the layering technique in the Carrara Interaction manual.
- The post-anatomical pressed restoration can be perfectly fired and finished using Carrara Interaction.

Tp:

With the application of large ceramic volumes it is advised to first apply a thin layer of Antagon & Carrara (Spray) Glaze on the Carrara Volumia sub-structure. (The alternative wax in combination with Glaze powder.) The adding ceramic material can be applied immediately on this glaze "wash-layer" and then fired according to the Carrara Interaction denture program.

Colour-schedule for Carback-technique only

Indication	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Liner	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Pellet	T01	A2	A3	A3.5	A4	T04	B2	B3	B4	T02	T05	T03	C4	T06	D3	D4
Action-dentist	1&2	3&4	1&2	3&4	1&2	3&4	1&2	3&4	1&2	3&4	1&2	3&4	1&2	3&4	1&2	3&4
X-tra Incisals	Red - Blue - Orange - Grey - Clear															
Incisal	5&8	5&8	5&9	5&9	60	5&7	5&9	5&9	60	5&9	60	60	60	5&9	60	5&9

10. "Correktion firing"

- Carrara Interaction Correction (Art. no. 484455081) can be used for correcting firing.
- The firing instructions can be found in the Carrara Interaction manual.
- It is also possible to use Carrara Interaction Margin Correction (Art. no. 482125018 - ...020) for margin corrections.

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