
# Carrara Press

Persovens  
Press furnaces

Fours de pressée  
Pressöfen

English

Deutsch

Elephant

dental health products

step-by-step

# Index

- I** Preparation
- II** Investing (Silicon Muffle Ring)
- III** Investing (Muffle Paper)
- IV** Pre-heating (rapid process)
- V** Pre-heating (slow process)
- VI** Pressing
- VII** Finishing
- VIII** Carrara Press Core  
(preparation for baking with Carrara Vincent)
- IX** Carrara Press Inlay  
(finishing pressed objects)
- X** Appendix

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## Carrara Press "Step by step"

Carrara Press "Step by step"  
The full ceramic addition to the  
Carrara Premium System

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**The Carrara System:**  
perfect results through quality products

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step-by-step

# Preparation

- Prepare a working model using a (super-hard) die-stone.
- Apply at least two or three coats of die spacer on to the prepared die.
- Leave at least 1 mm of the margin clear of die spacer.
- The optimum final coating thickness is approx. 30 µm.

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- To prepare the required wax model, use a fully residue-free combustible modelling wax (see appendix).

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
**Dimensioning:****Carrara Press Core (caps)**


- The minimum thickness after processing the Carrara Press Core caps is 0.8 mm.

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**Carrara Press Inlay (objects)**

- The minimum thickness of Carrara Press Inlay objects should be 0.8 mm measured at the thinnest point (see dentists' preparation guidelines).

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- Seal the opening at the tip of the Former with wax.

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- Choose the appropriate Former.

|              |                 |                    |
|--------------|-----------------|--------------------|
| 100 g former | 1 object        | 0,2 - 0,6 g of wax |
| 200 g former | 2-4 objects     | 0,2 - 1,2 g of wax |
|              | Max. 3 molars   |                    |
|              | Max. 4 incisors |                    |

- Use a  $\varnothing$  3.5 mm wax channel with a length of 8 mm.
- Clean the Former thoroughly and remove any traces of wax.

- Carefully apply the press channel (see illustration 7).
- Ensure that all transitions blend.
- Avoid sharp edges.

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
8

8

### Weight of wax, including Carrara Press Channel

- Weigh the wax objects, including the press channel.

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- Check the correct position of the object.  
(see illustration 7).

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# Investing

## Silicon Muffle Ring

Concerns the standard Carrara Press Investment, article number 232768040.

For instructions on the use of Carrara Press Speed Investment, article number 232772040, please see the user's manual supplied with the investment.

- Spray the waxed-up objects with an alcohol-free surface-tension-reducing agent.

- Afterwards, carefully blow-dry the objects.

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- Treat the Silicon Muffle Ring with silicon spray (oilfree), which facilitates removal of the ring from the investment material.

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- Place the Silicon Muffle Ring on the muffle base.

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- Mix the Carrara Press Investment under vacuum for 60 seconds (100 g powder with 24 ml diluted liquid).
- The recommended dilution ratio of Carrara Press Investment Liquid to distilled water is 50-50% up to 60-40%

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- Immediately after filling, attach the Muffle Topshaper with a twisting movement.
- After filling, allow the CP Investment to set for 15 minutes.
- If the setting time is exceeded, the muffle can no longer be heated quickly. In that case, heat the muffle step by step (slow process).

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- Remove the Muffle Topshaper and Silicon Muffle Ring.

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- Remove any excess investment material using a scalpel.

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- Shortly before removing the muffle base, slightly roughen the underside of the muffle with sandpaper (120 grit). Make certain that the underside remains even.

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- Remove the Muffle Base from the muffle using a slight twisting movement.
- Remove the upright investment edge.
- If necessary remove any investment material residue with compressed air.

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# Investing

## Muffle Paper

Concerns the standard Carrara Press Investment, article number 232768040.

For instructions on the use of Carrara Press Speed Investment, article number 232772040, please see the user's manual supplied with the investment.

- Treat the wax objects with a surface tension reducing agent.

- Afterwards, carefully blow-dry the objects.

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- Attach the dark blue edge of the ringformer paper on the desired marking (for 100 g applications, trim the excess paper along the dotted line).

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- Position the ringformer paper on the Muffle Base and then attach the stabilisation ring.

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- Mix the Carrara Press Investment under vacuum for 60 seconds (100 g powder with 24 ml diluted liquid).
- The recommended dilution ratio of Carrara Press Investment Liquid to distilled water is 50-50% up to 60-40%.

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- Fill the investment into the muffle on a model vibrator.

- Remove the stabilisation ring and attach the Muffle Topshaper using a twisting movement.
- After filling, allow the Carrara Press Investment to set for 15 minutes.
- If the setting time is exceeded, the muffle can no longer be heated quickly. In that case, heat the muffle step by step (slow process).
- Remove the Muffle Topshaper and the paper.

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- Remove the seam that is formed in the muffle wall using a scalpel.
- Shortly before removing the muffle base, slightly roughen the bottom of the muffle with sandpaper (120 grit). Make certain that the bottom remains even.

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- Remove the Muffle Base from the muffle using a slight twisting movement.
- Remove the upright investment edge.
- If necessary, remove any investment material residue with compressed air.

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# Pre-heating

Rapid process

(Only possible if the setting time of 15 minutes is not exceeded.)

### Direct heating (rapid process):

- 1. Place the muffle in a pre-heated furnace at a temperature of 600 °C.
- 2. The muffle should stand in the furnace at a temperature of 600 °C for 10 minutes.
- 3. Allow the temperature to rise to 950 °C over a period of 15 minutes (10 °C/min.).
- 4. Leave the muffle to stand at 950 °C for at least 45 minutes.
- 5. The maximum time in the pre-heating furnace must not exceed 3 hours.

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

Slow process

## Phased heating (slow process):

- 1. Place the muffle in a cold furnace.
- 2. Allow the temperature of the furnace to rise at a rate of 5 °C/min. until it reaches 250 °C.
- 3. Maintain the temperature at 250 °C over a period of 30 minutes.
- 4. Then increase the temperature at a rate of 5 - 10 °C/min. until it reaches 950 °C.
- 5. Leave the muffle to stand at 950 °C for at least 45 minutes.

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# Pressing

- Take the pre-heated muffle out of the furnace.
- After about 60 seconds, place the required quantity of Carrara Press Pellets, one by one in the muffle duct (minimum 2, maximum 4).
- Make sure you have the right type and shade (see table 6 and 7).
- Carefully place the pellets in the press channel, making certain that they do not get stuck.

**Weight of wax, including Carrara Press Channel**  
Always use at least 2 Carrara Press Pellets.

- Up to a max. of 0.6 gr : wax: 2 Carrara Press Pellets
- Up to a max. of 0.9 gr : wax: 3 Carrara Press Pellets
- Up to a max. of 1.2 gr : wax: 4 Carrara Press Pellets

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- Insert the Carrara Press Plunger (37 mm) with the smooth side facing down.

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- Place the muffle in the furnace immediately and press according to the appropriate furnace programme (see pressing table).

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Empress and EP 500 are registered trademarks of Ivoclar Dental.

- Remove the muffle from the furnace, using the muffle tongs.
- Let the muffle cool down completely to room temperature.

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- The caps for manufacturing crowns are pressed with Carrara Press Core Pellets.

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- Inlays, onlays and veneers are pressed using Carrara Press Inlay Pellets.

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## Finishing

- Devest the castings by sandblasting with glass beads (50 µm).
- Maximum pressure 2.5 bar!

- Carefully remove the castings from the press channels using a diamond disc (not diaflex).
- To avoid localised overheating of the ceramic, make sure that neither the working pressure of the disc nor the rotational speed (max.10,000 r.p.m) is too high.

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- Finish the ceramic object and, if necessary, shape it using sharp and (very) fine-grained diamond burs. (max. 15,000 rpm).
- Note the prescribed wall thickness! (see chapter I, point 3 and 4, Dentists' preparation guidelines).

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- Remove the die spacer and fit the object carefully onto the die.
- Identify any defects using contact fluid.
- Remove any high spots with a fine diamond bur (max. 10,000 rpm).

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
dental health products

# Carrara Press Core

## Preparation for baking with Carrara Vincent

- The Carrara Press Core objects are baked with Carrara Vincent Ceramic.
- First apply a layer of conditioner to the Carrara Press Core caps.

- Mix the Carrara Press Conditioner with the Elephant Stain Liquid until you have a thin consistency.
- Apply a very thin layer ("Wash firing") of the mixture on to the cap.

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- Place anteriors, premolars and molars on a firing pillow tray.

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- Fire the conditioner in accordance with the table presented above, without a cooling phase. Mix the conditioner with Stain Liquid.
- Perform each subsequent firing programme in accordance with the usual Carrara Vincent firing procedure (including OD!), without a cooling phase.
- Glaze the objects (Universal Glaze!) in the same way as Carrara Vincent objects.

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### Possible shape corrections:

- The Carrara Press Inlay objects cannot be baked with Carrara Vincent.
- Only use Carrara Correction material to alter the shape of the Carrara Press Inlay objects.

This complementary material is used for:

- aesthetic/functional occlusal/incisal adjustment;
- optimising the surface structure to increase aesthetic value;
- corrections, such as contact points, etc.

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**Working procedure:**

- Shape your wax-up in the desired (mamelon) form and press the object as usual.
- Or: Reduce the previously pressed object in the desired areas with diamond finishing burs.
- Note: Be aware of overheating during grinding (max. 15,000 rpm).
- Rinse / clean the surface thoroughly before applying the desired material(s).
- Warning: Do not use conditioner.
- Mix the powder and carving liquid to a regular consistency.
- Apply the material at the areas needed and allow it to dry sufficiently.
- During firing, the objects need to be supported by a thin layer of ceramic firing pillow.

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## Glaze firing:

- Carrara Press Inlay objects need to be glazed with Carrara Press Glaze
- Depending on the required lustre, one or two coatings of Carrara Press Glaze will be sufficient.
- Mix the powder with Stain Liquid.
- Apply the Carrara Press Glaze material (and Stains, if required) and fire according to the firing chart (see below).

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## Appendix

Recommended wax types for application with Carrara Press:

- Aldente - pro-mod VKS, grey
- pro-mod modelling wax for K+B, gray, normal
- pro-mod modelleerwas voor K+B, green, hard
- pro-mod modelling wax for K+B, green, regular
- Benzer Dental - Chromowax Nr.0 - Nr.9
- DT - Finocrown, grey
- Finowax dipping wax, green
- Renfert - Crowax, blue for K+B
- Pico, beige, for K+B
- Schuler - Modelling wax, blue
- Yeti - VKS modelling wax, beige-transparent
- VKS dipping wax, bordeaux red
- VKS modelling wax, grey-transparent
- Dentaorium - Spezial-Gusswachs by Gründler, red
- Dipping wax, orange
- Special modelling wax for press ceramic, beige

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