### Working procedure for BIOGOLD INTERACTION

- 1. Indication.** Rich yellow alloy for metal-ceramic, inlays, crowns, bridges and telescopic technique.
- 2. Contraindication.** If patients are known to be allergic to any of the components, the material should not be applied.
- 3. Sprueing.** According to the indirect method: 3 mm  $\varnothing$  for crowns, 3,5 mm  $\varnothing$  for pontics, runner bar 4-4,5 mm  $\varnothing$ .
- 4. Investing.** Invest in a reliable phosphate-bonded investment, such as Carrara Universal DL Investment.
- 5. Burn out.** Preheat at 800°C/1470°F for 30-120 min. depending on the size of the muffle ring.
- 6. Casting.** Casting temperature: 1235°C/2255°F in a graphite crucible; or with the flame in a ceramic crucible with for example Elephant Proflux. Add at least 1/3 new metal at each casting. Clean the metal well before re-use, for example by means of sandblasting and steam-cleaning.
- 7. Pickling.** 5-10 Minutes in warm Pickling Agent.
- 8. Finishing.** Use carbide burs and/or ceramic bonded stones. The final trimming must be carried out in one direction only by means of fine crossteethed carbide burs! Sandblast with 50-125  $\mu$ m aluminum oxide (Alublast), (pressure max. 2 bar).
- 9. Cleaning.** With distilled water ultrasonic or by means of a steam cleaner.
- 10a. Oxidation.** 3 Minutes at 950°C/1740°F without vacuum. Do not exceed a heat-rate of 55°C/131°F in order to avoid an override of the maximum temperature.
- 10b. Oxidreduction.** Sandblast with aluminum oxide 125  $\mu$ m (pressure max. 2 bar) or pickle for max. 1 min. in clean pickling agent. Afterwards clean again (see point 9).
- 11. Bake.** With a standard fusing and expandic ceramic, such as Antagon Interaction. Baking temperature may not exceed 950°C/1740°F. If the porcelain manufacturer prescribes a baking temperature higher than 950°C/1740°F, the baking time at 950°C/1740°F should be extended by one minute for every additional 10°C/50°F.
- 12. Pre-soldering.** Orion 1030°C/1885°F yellow PdF with a flux for example Elephant Paste Flux.
- 13. Post-soldering.** Elephant II PdF 800°C/1470°F or III PdF 750°C/1380°F with a flux for example Elephant Paste Flux.
- 14. Soft annealing (for inlays HV5-140).** 10 Minutes at 900°C/1652°F, then quench to casting object in water.

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