


User instructions



Available in the following colours.
margin 1A2, 3A4, 1B2, 3B4, 1C2, 3C4 & Booster
margin correction Bright, Medium & Dark

Shoulder preparation instructions


These instructions apply to Carrara Interaction ceramic shoulder preparations. For Standard Carrara Interaction preparations, please refer to the Carrara Interaction user instructions.



1 Sealing the die

1 Sealing the die:

- Seal the die by applying a die hardner agent and allow to dry well.




2 Separating the die

2 Separating the die:

- Apply a ceramic separation agent to the die and dry with compressed air.

| Colour | A1 | A2 | A3 | A3,5 | A4 | B1 | B2 | B3 | B4 | C1 | C2 | C3 | C4 | D2 | D3 | D4 |
|--------------------------|-----------------------|----|-----|-----------------------------------|---------|-----|----|-----|-----------------------------------|-----|----|-----|-----------------------------------|-----|-----|-----|
| <i>margin</i> | 1A2 | | 3A4 | 50% 3A4 & 50% Booster | Booster | 1B2 | | 3B4 | 50% 3B4 & 50% Booster | 1C2 | | 3C4 | 50% 3C4 & 50% Booster | 1A2 | 3A4 | 3C4 |
| <i>margin correction</i> | Bright, Medium & Dark | | | | | | | | | | | | | | | |

A close-up photograph showing a ceramic shoulder being prepared. A brush is applying a white, cream-like material to the shoulder area of a light-colored ceramic piece.

3a Applying *margin*

3 Applying *margin* for the first shoulder:

- Use **margin liquid** to mix the **margin** in question (see the colour schedule on the previous page) into a cream-like consistency.
- Apply the **margin** material in the usual manner. Before proceeding, condense the material and dry it with a hairdryer, if required.

A close-up photograph showing a ceramic shoulder being dried. A hairdryer is directed at the shoulder area, which is covered in a white, cream-like material.

3b Dry with hairdryer

A close-up photograph showing a ceramic shoulder being removed from a die. The shoulder is covered in a white, cream-like material.

4 Removing the object

4 Removing the object:

- Remove the object from the die.
- Fire the ceramic in accordance with the firing chart.

4 Firing chart for the first shoulder preparation firing:

| Drying | Starting temp | Vacuum- starting temp. | Rate of heat increase |
|-------------|----------------|------------------------|-----------------------|
| 8–10 Min. | 450°C | 450°C | 55–80°C/Min. |
| Final temp. | Retention time | Surface appearance | |
| 860°C | 3 Min.# | Matte finish | |

2 min. with vacuum and 1 min. without vacuum

A close-up photograph showing a ceramic shoulder after firing. The shoulder is covered in a white, cream-like material.

5 Shoulder after firing


5 Applying *margin* for the second shoulder preparation firing:

- Use **margin liquid** to mix the **margin** into a cream-like consistency.
- Apply the **margin** material in the usual manner. Before proceeding, condense the material and dry it with a hairdryer, if required.
- Remove the object from the die.
- Fire the ceramic in accordance with the firing chart.

5 Firing chart for the second shoulder preparation firing:

| Drying | Starting temp | Vacuum- starting temp. | Rate of heat increase |
|-------------|----------------|------------------------|-----------------------|
| 8–10 Min. | 450°C | 450°C | 55–80°C/Min. |
| Final temp. | Retention time | Surface appearance | |
| 860°C | 2 Min.# | Matte finis | |

1 min. with vacuum and 1 min. without vacuum



6 Applying margin correction

6 Processing:

Correcting the ceramic shoulder preparation using margin correction:

- Use **margin liquid** to mix **margin correction** into a cream-like consistency.
- Apply **margin correction** to the location to be corrected and condense it. Before proceeding, dry with a hairdryer.
- Remove the object from the die.
- Fire the ceramic in accordance with the firing chart.
- After the firing process, the correction layer should be polished mechanically, using silicone polishing instruments and diamond paste.

Pictures: ZTM Jan Schünemann (Bielefeld, Duitsland)

6 Firing chart for margin correction (Bright, Medium, Dark)


| Drying | Starting temp | Vacuum- starting temp. | Rate of heat increase |
|-------------|-----------------------|------------------------|-----------------------|
| 6 Min. | 450°C | 450°C | 55–80°C/Min. |
| Final temp. | Retention time | Surface appearance | |
| 795°C | 2 Min. without vacuum | Silky gloss | |

After the firing process, the correction layer should be polished mechanically.

Carrara Interaction: One ceramic for all Carrara System applications.


Ceramic



 Carrara Interaction


Universal alloys



 Carrara PdF
High gold content

 Cera F
Reduced gold content

 Cera R plus

 Cera E
Palladium-silver

Press ceramic



 Carrara Press
Full ceramic

One ceramic for both alloys and press ceramics: Carrara Interaction

Elephant

We support your success.

Elephant Dental B.V.
Verlengde Lageweg 10
1628 PM Hoor, The Netherlands
Tel. +31 229 25 90 00
Fax +31 229 25 90 99
E-mail export@elephant.nl
Internet www.elephant-dental.com

dental health products



0344

11-dec-2003