Carrara Universal DL Investment

gießen, casting, gieten, moulage, lanciare, lanzar

60 - 150 g

Elephant

dental health products

CUDL 1003

Inhaltsangabe
Index
Inhoudsopgave

Protocole
Indice
Indice

D

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FR

GB

3	Product profile
4	Physical properties
5	Remarks
6	Investing
11	Pre-heating
16	Finishing

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Product profile:

Carrara Universal DL Investment is a universal phosphate- bound, graphite- free dustless investment for the casting of dental precious alloys and pressing ceramics. The investment powder comes packed in sachets of 60, 100, 150 g or bucket of 6 kg. The associated Carrara Universal DL Investment Liquid is a concentrated mixing liquid, which, when diluted with distilled water, guarantees perfect workability and fit.

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GB

Physical properties:

Mixing ratio:	60 g/14,4 ml - 150 g/36 ml
Mixing time:	60-120 sec. under vacuum
Processing time:	7 Min. at 21°C/70°F
Max. expansion during setting:	ca. 2%
Max. thermal expansion:	1,1% (at 500°C/932°F)
Compressive strength:	6,7 MPa

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

Remarks:

1. The working time of Carrara Universal DL Investment can be increased by keeping the liquid in the refrigerator. Close the bottle after use. When frozen or evaporated, the liquid will crystallize on the bottom, and the liquid can no longer be used.
2. STORAGE: Store powder in a dry and cool place (temp.: 8-25°C/46-77°F).
3. Use the mixing bowl for phosphate-bonded investments only.
4. Carrara Universal DL Investment contains silica. This release is reduced to a minimum because of the patented treatment. Avoid inhalation of dust.
5. Ensure adequate extraction and ventilation during the burning-out process.

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GB

• Follow the manufacturer's instructions to apply the casting sprues. We recommend the indirect method.

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GB

- Line the casting rings with moistened Olympos-Liner (asbestos free).

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GB

- Mix powder and liquid with a spatula for 10-15 sec. The recommended dilution of Carrara Universal DL Liquid to distilled water is mentioned on page 20.

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GB

- Mix in a mechanical vacuum mixer for 60-120 sec. (at 360 rpm).

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GB

- Scrape away the skin on top of the casting ring and place it in a cold burnout furnace.

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
- Scrape away the skin on top of the casting ring and place the muffle in a burnout furnace at a preheating temperature of 600-950°C/1112-1742°F.

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
GB

- Raise temperature gradually (5-10°C/41-50°F per min.) to maximum pre-heat temperature (950°C/1742°F), depending on the alloy to be cast. Hold end temperature for 30-90 minutes before casting, depending on ring size.

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
GB

- Raise temperature to maximum pre-heat temperature (950°C/1742°F), depending on the alloy to be cast. Hold end temperature for 30-90 minutes before casting, depending on ring size.

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GB • Once the muffle has cooled down completely, place it in water or hold it under running water to avoid dust generation during devesting.

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GB

- The investment can be separated from the muffle ring more easily using a devesting press.

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For the above mentioned alloys, we recommend to use the directives for the following liquid-concentration:

50 – 55 % Carrara Universal DL Investment Liquid + 50 – 45 % distilled water or 18 – 19,8 ml Carrara Universal DL Investment Liquid + 18 – 16,2 ml distilled water.

45 – 50 % concentrate	Inlays, castingchannels
50 % concentrate	Crowns, bridges, conuscrowns (primary parts) and onlays
40 % concentrate	Telescopic crowns (secondary parts)

Elephant alloys group II

Cera Aurum	Minerva 3S
Orion WX	Minerva 4PdF
Orion E	Minerva 4CF

For these alloys we recommend to use the directives for the following liquid-concentration:

55 – 60 % Carrara Universal DL Investment Liquid + 45 – 40 % distilled water or 19,8 – 21,6 ml Carrara Universal DL Investment Liquid + 16,2 – 14,4 ml distilled water.

50 – 55 % concentrate	Inlays, castingchannels
55 % concentrate	Crowns, bridges, conuscrowns (primary parts) and onlays
45 % concentrate	Telescopic crowns (secondary parts)

Elephant alloys group III

Orion Vesta	Pallacon
Orion Isis Plus	Pallium 3C
Orion Argos	Pallium 3
Orion Virgo P	Cera E

When a ringless casting technique is used, we recommend reducing the liquidconcentrates mentioned above by 10-20 %.

Quick process:
starting temperature 600°C/1112°F.

For the alloys of group III, we recommend to use the directives for the following liquid-concentration:

65 – 70 % Carrara Universal DL Investment Liquid + 35 – 30 % distilled water or 23,4 – 25,2 ml Carrara Universal DL Investment Liquid + 12,6 – 10,8 ml distilled water.

60 – 65 % concentrate	Inlays, castingchannels
65 % concentrate	Crowns, bridges, conuscrowns (primary parts) and onlays
55 % concentrate	Telescopic crowns (secondary parts)

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