# PERFLEX<sup>®</sup>

## PERFLEX DUPLICATING MATERIAL Step by Step Instructions

PERFLEX is a reversible hydrocolloid duplicating material which has been designed to produce refractory models (V.R. Investment). It may be used to produce molds for the fabrication of partial and full dentures utilizing the fluid resin casting technique.

### Melting Procedures:

Remove cover from the container and remove as much material as required. Cover and seal any remaining portion in the original container.

Chop material into small pieces; heat and carry out the melt cycle in a double boiler steam pot. Avoid the use of aluminum, copper and iron bearing alloyed containers. If automatic melt equipment is employed, follow manufacturer's instructions. However, the melt temperature should be set at 97°C (207°F) for no more than ONE (1) HOUR. Remove from source of heat and allow the melt to cool 55°C (131°F) while agitating. Maintain the melted material at a POURING TEMPERATURE OF 55°C (131°F). This precaution will extend the life of the material.

IF A HOTPLATE IS USED, exercise care to avoid overheating and/or scorching. Stir while heating until all the particles are melted and completely in solution. Heat the melt to 97°C (207°F) for no more than ONE (1) HOUR. Remove from source of heat and allow the melt to cool 55°C (131°F) while agitating. Maintain the melted material at a POURING TEMPERATURE OF 55°C (131°F). This precaution will extend the life of the material.

### Usage:

Fill flasks containing properly soaked models with PERFLEX Duplicating Material and chill in cold water bath 10°C (50°F ± 10°F) for 45 minutes. Chill longer if water temperature is warmer than 15°C (60°F).

### Refractory Model Duplication:

Remove model and gently blow dry with air jet and place mold (open side down) on bench to retard drying and shrinkage. Follow manufacturer's instructions for mixing and setting times of the refractory slurry.

### Stone Model Duplication:

The procedure for mold pouring and preparation (see Refractory Model Duplication) is observed for Stone Model Duplication. After blow-drying the mold, fill the mold cavity with PERFLEX Mold Soak ReOrder# NO26535, and allow to set on bench for fifteen (15) minutes. Return the wash to its container and blow dry the mold with an air jet. Mix stone as per manufacturer's instructions. Slower setting stones usually require a longer period of spatulation to achieve a faster initial set.

Fast setting stones usually require a minimum of 45 minutes dwell before removal; whereas the slower setting stones will require a 60 minute mold dwell time. At the proper time, remove mold with duplication from the flask and cut mold away from duplicate model. Trim the model and rinse with water.

### Fluid Resin Casting:

Properly soak model with wax pattern and position in a flask in desired manner. Fill the flask with molten material. Chill in a cold water bath 10°C (50°F ± 10°F) for forty (40) minutes. Position flask in the water bath on slats or ribs to allow proper water circulation. After chilling, disassemble the flask and remove the mold. Silt mold in three of four places (being careful not to cut down on the waxed area), and remove model. Cut desired sprues and venting holes and rinse mold in rapidly running cold water. Blow dry with an air jet and cover mold not being worked on with a damp towel. Follow resin manufacturer's technique for casting of dentures.

### Notes:

After molds have been used for either model duplication or fluid resin casting, they should be rinsed free of any stone or refractory debris, and stored in a CLOSED, AIRTIGHT container until the next remelt cycle. Any liquid that separates from the used molds during storage should NOT be discarded but added to the mold pieces during the next remelt cycle.

After several remelts, the duplicating material may thicken and pour slowly at the pouring temperature of 55°C (131°F). The resulting molds may appear dry and somewhat brittle. To eliminate this problem, add 150 ml of water for EACH bucket of PERFLEX being melted at the beginning of the remelt cycle.

Grainy PERFLEX may be obtained if the material is not heated to 97°C (207°F) and maintained for a sufficient period of time to dissolve all the particles.

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## MATERIAL SAFETY DATA SHEET

FORM NO: 314 SECTION 1 - GENERAL INFORMATION

1.1 Product Trade Name.....: Perflex Duplicating Material  
1.2 Part (Item) Number.....: N021055  
1.3 Chemical Name.....: NA  
1.4 Synonyms.....: Hydrocolloid  
1.5 Division Name.....: Dentsply Prosthetics  
Address.....: 570 W. College Ave.  
City, State, Zip.....: York, PA 17405-0872  
1.6 Phone No. for Info.....: 717-845-7511  
"FOR CHEMICAL EMERGENCY"  
SPILL, LEAK FIRE, EXPOSURE OR ACCIDENT  
CALL CHEMTREC - DAY OR NIGHT  
1.7 Emergency Phone No.....: 1-800-424-9300 1.12 National Stock No: NA  
CHEMTREC  
1.8 Date.....: 2/22/05 1.13 Specification No: NA  
1.9 Dentsply Contact.....: Susan M. Hanson 1.14 Chemical Family: NA  
1.10 Signature.....: *Susan M. Hanson* 1.15 Type/Grade/Class: NA  
1.11 Formula.....: NA 1.16 NRC/EPA Number: NA

## SECTION II - HAZARDOUS INGREDIENTS

2.1 CAS Number Chemical Name (Ingredients) % TLV OSHA PEL  
NONE

## SECTION III - PHYSICAL PROPERTIES

3.1 Boiling Point: NA 3.11 Specific Gravity: 1.10 @ 25°C  
3.2 Vapor Pressure: NA 3.12 pH.....: NA  
3.3 Vapor Density: NA 3.13 Evaporation Rate: Slight  
3.4 Critical Temp: NA 3.14 Viscosity: NA  
3.5 Decomposition Temp: NA 3.15 % Volatile by Vol: NA  
3.6 Melting/Freezing Pt: 98°C 3.16 Magnetism: NA  
3.7 Solubility in Water: Slight 3.17 Autoignition Temp: ND  
3.8 Critical Pressure: NA 3.18 Corrosion Rate: NA  
3.9 Perm. Exp. Limit: ND 3.19 Molecular Weight: NA  
3.10 Appearance and Odor: Translucent brown semi-solid gel.

## SECTION IV - FIRE AND EXPLOSION HAZARD DATA

4.1 Flash Point (Method Used).....: Non-flammable.  
4.2 Flammable (Explosive) Limits :  
LEL: NA UEL: NA  
4.3 Extinguishing Media.....: Non-combustible, use extinguishing media  
appropriate to the surrounding fire conditions.  
4.4 Special Fire Fighting Procedures and Equipment: None.  
4.5 Unusual Fire and Explosion Hazards.....:None.

## SECTION V - HEALTH HAZARD DATA

5.1 Effects of Acute and Chronic over Exposure: Eye contact may result in slight irritation.  
5.2 Emergency and First Aid Procedures: In case of contact with eyes, flush with large quantities of water for 15 minutes. Contact a physician.

## SECTION VI - REACTIVITY DATA

6.1 Unstable: Stable: X  
6.2 Conditions to avoid (Stability): None.  
6.3 Incompatibility (Materials to avoid): None.  
6.4 Hazardous Decomposition Products: None  
6.5 Hazardous Polymerization: May Occur: Will Not Occur: X  
6.6 Conditions to avoid polymerization: None.

## SECTION VII - SPECIAL PRECAUTIONS AND SPILL OR LEAK PROCEDURES

7.1 Handling and Storage Precautions: Avoid getting material in eyes.  
7.2 Other Precautions: None.  
7.3 Steps to be taken in case material is released or spilled: If material is spilled while molten, allow it to cool and solidify. Place solid material in appropriate container.  
7.4 Waste Disposal Procedures: Dispose of material in approved landfill. Check current local, state and federal regulations before disposal.  
7.5 Neutralizing Agent: NA

## SECTION VIII - OCCUPATIONAL PROTECTION MEASURES

8.1 Respiratory Protection: None.  
8.2 Ventilation: None.  
8.3 Protective Gloves: None.  
8.4 Eye Protection: Not required.  
8.5 Other Personal Protective Equipment: None.

## SECTION IX - TRANSPORTATION

9.1 Applicable Regulations: IMOC: IATA: 49 CFR:  
Military Air (AFR 71-4): Tariff 6D:  
9.2 Shipping Name.....: NA  
9.3 I.D. Number.....: NA 9.9 Unit Container.....: NA  
9.4 Report Quantity.....: NA 9.10 DOT SPSC Container.....: NA  
9.5 Hazard Class.....: NA 9.11 DOT Exempt/DOD/CCN.....: NA  
9.6 Labels.....: NA 9.12 Limited Quantity.....: NA  
9.7 Aerosol Propellant(s).....: NA  
9.8 Net Explosive Weight.....: NA

**DENTSPLY**  
AUSTENAL

DENTSPLY International Inc.  
DENTSPLY Austenal  
York, PA 17405-0872  
1-800-621-0381  
www.austenal.dentsply.com  
www.dentsply.com

Perflex® is a Registered Trademark of DENTSPLY Austenal.  
Made & Printed in USA  
NPK0814 Rev. A (9/05)