DENTSPLY
AUSTENAL

Vi-COMP[®] Investment

DIRECTIONS FOR USE

Introduction:

A high expansion/high heat investment especially formulated for use when casting Vi-Comp Alloy. Selective dilution of special liquid allows the user to control expansion.

Spruing

Indirect spruing is recommended for multiple units; for single units, a direct sprue is used. Use a 10 gauge sprue for copings and a 6 gauge sprue for bulky crowns and pontics. Sprues should not be less than 6.5mm and not more than 9mm in length. Copings should be waxed not less than 0.5mm thick.

Investing

Because metal rings constrain setting expansion, the use of paper rings (ReOrder #N123800) is highly recommended.

The expansion of Vi-Comp Investment is controlled by dilution of the special liquid. Maximum expansion is obtained by using full strength liquid and paper rings. If less expansion is needed, dilute the special liquid with water, but do not change the liquid to powder ratio.

Center the paper ring around the pattern and wax in place. If a wetting agent (debubblizer) is used, be sure to shake any excess out of the ring before investing.

- Use 10ml of liquid (full strength) for each 60g envelope
- Incorporate powder in liquid and then vacuum mix at slow speed for 30 seconds.
- If hand mixing, spatulate for 60 seconds
- Slowly pour mix in the ring while vibrating.
- Allow the mold to bench set for one (1) hour. Do not store overnight in humidior as this may result in cracking.

Burnout

Burnout should be done the same day if possible. (If burnout is not done the same day, soak the mold in room temperature water for 5 minutes prior to burnout. After soaking allow the mold to dry for 10 or 15 minutes before placing in a cold furnace.)

Place the molds in a cold furnace and raise the temperature to 1800°F (982°C) in less than 1¹/₂ hours. Allow a one hour heat soak at 1800°F (982°C).

Casting

Best results are achieved by minimizing the amount of time between removal of the mold from the burnout furnace and casting the metal. Allow molds to cool until they are warm to the touch before breaking out castings.

Caution

Special liquid is stable at room temperature, however, if the liquid freezes, it is unusable.

If any container you receive has any significant amounts of sediment on the bottom, do not use it.

During cold weather, liquid may freeze in transit, so please fill your winter requirements early.

Use separate mixing bowls for phosphate and gypsum investments.

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Vi-Comp® is a registered trademark of DENTSPLY Austenal.
©2005 DENTSPLY International Inc. All rights reserved. Printed in U.S.A. A1003-C Rev. A (8/05)
DENTSPLY International Inc.