Working procedure for ORION ISIS

- 1. Indication.** White dental alloy for metal-ceramic, inlays, crowns, bridges and telescopic technique.
- 2. Contraindication.** If patients are known to be allergic to any of the components, the material should not be applied.
- 3. Sprueing.** According to the indirect method: 3 mm \varnothing for crowns, 3,5 mm \varnothing for pontics, runner bar 4-4,5 mm \varnothing .
- 4. Investing.** Invest in a reliable phosphate-bonded investment, such as Carrara Universal DL Investment.
- 5. Burn out.** Preheat at 900°C/1650°F for 30-120 min. depending on the size of the muffle ring.
- 6. Casting.** Casting temperature is 1400°C/2550°F in a ceramic crucible with Elephant Proflux. Use at least 1/3 new metal at each casting. Clean the metal well before re-use, for example by means of sandblasting and steam-cleaning.
- 7. Finishing.** Use carbide burs and/or ceramic bonded stones. The final trimming must be carried out in one direction only by means of fine crossteethed carbide burs! Sandblast with 50-125 micron aluminium oxide (Alublast), (pressure max. 3 bar).
- 8. Cleaning.** Ultrasonic or steam cleaner.
- 9a. Oxidation.** 5 minutes at 950°C/1740°F, in atmosphere.
- 9a. Oxireduction.** Blast lightly with aluminiumoxide 50-125 micron (Alubast) until oxide is just removed. Afterwards clean again (see point 8).
- 10. Bake.** With a standard fusing ceramic, such as Antagon Interaction.
- 11. Pre-soldering.** Orion 1090°C/1990°F white with Elephant Paste Flux.
- 12. Post-soldering.** Elephant PP/820°C/1470°F with Elephant Paste Flux.

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