### Working procedure for ORION UX PLUS

- 1. Indication.** White dental alloy for metal-ceramic, inlays, crowns, bridges and telescopic technique.
- 2. Contraindication.** If patients are known to be allergic to any of the components, the material should not be applied.
- 3. Sprueing.** According to the indirect method: 3 mm  $\varnothing$  for crowns, 3,5 mm  $\varnothing$  for pontics, runner bar 4-4,5 mm  $\varnothing$ .
- 4. Investing.** Invest in a reliable phosphate-bonded investment, such as Carrara Universal DL Investment.
- 5. Burn out.** Preheat at 800°C/1470°F for 30-120 min. depending on the size of the muffle ring.
- 6. Casting.** Casting temperature is 1350°C/2460°F in graphite crucible; or with the flame in a ceramic crucible with Elephant Proflux. Use at least 1/3 new metal at each casting. Clean the metal well before re-use, for example by means of sandblasting and steam-cleaning.
- 7. Finishing.** Use carbide burs and/or ceramic bonded stones. The final trimming must be carried out in one direction only by means of fine cross-teethed carbide burs! Sandblast with 50-125 micron aluminium oxide (Alublast), (pressure max. 3 bar).
- 8. Cleaning.** With distilled water ultrasonic or by means of a steam cleaner.
- 9. Oxidation.** 5 minutes at 950°C/1740°F, without vacuum.
- 10. Bake.** With a standard fusing ceramic, such as Antagon Interaction.
- 11. Pre-soldering.** Orion 1060°C/1946°F (PdF) yellow or 1090°C/1994°F white with Elephant Paste Flux.
- 12. Post-soldering.** Elephant II (PdF) 800°C/1470°F or III (PdF) 750°C/1380°F with Elephant Paste Flux.

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