Dental ceramic system for over-pressing zircon oxide sub-structures for crowns and bridges.

Rx only
 ISO 6872

Elephant Dental B.V.
 Verlengde Lageweg 10
 1628 PM Hoorn, The Netherlands
 Tel.: +31 229 25 90 00
 Fax: +31 229 25 90 99
 E-mail: info@elephant.nl
 www.elephant-dental.com

Made in Holland

1. General information

1.1 Product description
 Sakura Volumia is a new system that offers an interesting and economical alternative to the dental ceramic layering technique. The Sakura Volumia press-technique is designed to press ceramic onto zircon oxide sub-structures using the 'full contour' technique.
 The combination of the colour-related *Linear* and the press pellets which transport harmonises between the dentin and the enamel enables the technician to make a reliable and simple reconstruction of the 16 classical 'V' colours. Sakura Volumia can also be used in the Carbock-technique.

1.2 Indication
 Dental ceramic for over-pressing zircon oxide sub-structures for crowns and bridges.
 Sakura Volumia can be used on zircon oxide sub-structures with a thermal expansion coefficient between 10.0 and 10.6 µm/m.K (25-500°C). Maximum wax weight (weight without sub-structure) of 1.4 g with a minimum strength of 0.8 mm (wax-up).

1.3 Contraindication
 The over-pressing in combination with:
 • Metal sub-structures
 • Aluminium oxide sub-structures
 • Glass infiltrated aluminium oxide sub-structures (Vita In-Ceram Zirconia)
 • All other, not mentioned sub-structures
 Bruxism and other parafunctions
 If patients are known to be allergic to one of the components, the material should not be applied

1.4 Reuse
 Sakura Volumia cannot be reused.

1.5 Precautions
 Consult the Material Safety Data Sheet (MSDS) for more information.

2. Preparation

2.1 Preparation
 To ensure a successful restoration using the Sakura Volumia system, please follow the general guidelines for preparation and layering thicknesses, as also described for full-ceramic restorations.

Preparation form	Circular preparation	Angle	Incisal/occlusal
Shoulder	1.5 mm (-/ 0.5 mm)	10°- 30°	1.5 mm - 2.0 mm
Chamfer	1 mm (-/ 0.5 mm)	10°- 30°	1.5 mm - 2.0 mm
Knife-edge	1 mm	10°- 30°	1.5 mm - 2.0 mm

3. Processing

3.1. Preparation of zircon oxide framework
 • The anatomical shape of the restoration should be considered when designing the sub-structure.
 • For the correct preparation of zirconium sub-structures, keep to the recommended process-steps in the user instructions from the supplier.
 • Afterwards apply a thin, masking layer of Sakura Volumia *Linear* and fire according the firing schedule below (see 3.2)
 • If required, a second coat can be applied (again afterwards).

3.2 Linear firing chart

Dry	Start temp.	Rate of heat increase	Final temp.	Vac.	Hold time
	400°C	60°C/min	920°C	1 min.	2 min.*

* 1 minute with vacuum, 1 minute without vacuum.
 Cooling: 0 minutes.

3.1. Wax-up
 • Before starting the wax-up, the restoration (incl. *Linear*) has to be weighed to determine the exact weight of the wax.
 • Prepare the required wax model 'full contour' by using a fully redwax-free combustible modelling wax suitable for press-ceramic system.
 • Keep in mind that the minimum thickness of the wax model should be as less than 0.8 mm.
 • With the Sakura Volumia technique, shape and occlusion of the wax-up can be perfectly modelled and checked in an articulator.

4.2 Layer thickness
 Layer thickness sub-structure min. 0.4 mm after finishing
 *To assure the colour aesthetic please keep to equally divided layer thickness among the bridge-units.

4.3 Spruing
 It is important to keep the following guidelines in mind when spruing the press-charset:
 • When spruing with only one press-charset, the diameter of the wax wire should be at least 2.5 mm Ø.
 • Do not exceed the length of max. 6 mm.
 • Large molars: attach a press charcoal of every coup (Ø 3 mm) to distribute the pressure evenly during pressing.
 • Preserve the occlusal contacts (IJJ buccal, U lingual).
 • The muffle former must be kept completely clean and free of wax or other residues.
 • Check the correct position of the object (see picture 1&2).
 • Avoid sharp edges (see picture 1 & 2).
 • Ensure that all transitions blend (see picture 1&2).
 • The spruing angle from the wax object to the muffle former must be 45° (see picture 1&2).
 • Determine the final wax weight, including the press-charset(s).

4.4 Wax weight chart including press channel

	1 pellet	2 pellets
Wax weight	up to 0.7 g	0.7 g - 1.4 g
Ring-size	100 g	200 g
Number of objects	Single crowns	Bridges or multiple crowns*

* Do not exceed the max. wax weight.

5.1 Investing
 For investing, use Carrara Universal Duettes Investment (art. no. 232978050) and follow the user instructions. (Recommended liquid-concentration: 50% Carrara Universal DI Liquid - 50% distilled water)

5.2 Burnout program (preheating)*

Start temp.	Rate of heat increase	Final temp.	Hold time
700°C	60°C/min.	940°C	30 min.

* Continuous maximum vacuum

* Please follow the manufacturer's instructions for your specific furnace.
 • Minimal burnout time for the Sakura Volumia system: 30 minutes.
 • Do not exceed the maximum burnout time of 2 hours!

*** Note:**
 • For each additional muffle, the minimal hold-time should be increased by 10 minutes.
 • The maximum hold-time should not exceed 2 hours!

6.1 Pressing
 • Take the muffle out of the furnace immediately after preheating.
 • Insert the required quantity of press pellets (see chart 4.4) and the Carrara Press Plunger (art. no. 004100580) with the smooth side facing down and place the muffle in the furnace.

Colour-schedule for Full-Contour only

Indication	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
<i>Linear</i>	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
<i>Transpa Dentin Pellet</i>	T01	T01	T02	T02	T03	T04	T04	T05	T05	T01	T02	T02	T03	T01	T06	T06
<i>Sakura Shade point</i>	1	2	2	2	2	4	5	4	4	5	6	7	7	3	8	9

*** Note:**
 • Please be aware that after inserting the first press pellet, it is advised to start by positioning the Carrara Press Plunger to ensure that the press pellet is correctly inserted.
 • Then insert the second press pellet and secure it with the Carrara Press Plunger, as before.

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1.4 Finishing
 • Carefully cut off the press channels with a diamond disc.
 • To avoid local overheating of the ceramic, make sure that neither the working pressure nor the rotational speed (max. 8.000 rpm) of the disc is too high.
 • The final finishing and shaping of the ceramic restoration can be done using sharp and fine-grained diamond coated discs (max. 15.000 rpm).
 • The entire surface of the restoration must be worked.

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