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Pre Solder Flux:bestelnummer

N200122, EU8200122

Fine Low Fusing Gold solder :bestelnummer

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CE
0086

R_Xonly **MSDS**

MSDS is available on the web or by request.

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Warning: This product contains substance(s) known to the State of California to cause cancer, or other reproductive harm. (Calif. Prop. 65)

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Vi-Comp® is a registered trademark of DENTSPLY Austenal.

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©2006 DENTSPLY International Inc.
All rights reserved.
Printed in U.S.A.
NPK1933 Rev. E (5/07)

USA **GB**

Indications
Vi-Comp II Alloy is a coball-based dental casting alloy, intended for metal ceramic or full metal dental restorations.

Contraindications

- Vi-Comp II Alloy is contraindicated for patients and users with a history of sensitivity to cobalt, chromium, tungsten, gallium, ruthenium and copper.

Warnings

- Prolonged exposure to alloy dust and/or fumes may lead to lung irritation, and/or pulmonary complications. Use appropriate engineering controls to limit exposure. For excessive inhalation of dust or fumes, seek medical advice.
- This alloy contains chromium. Some compounds of this element are potential carcinogens. Advise your health care provider of exposure to this element.
- Ingestion of gallium may cause bone marrow abnormalities with damage to blood forming tissues.
- This alloy contains elements that are known to the State of California to be carcinogenic.

Precautions

- When melting, grinding, or polishing alloys, use adequate ventilation, vacuum systems, protective eyewear, and protective masks and clothing.
- Check with the porcelain manufacturer to assure matched average thermal expansion of selected porcelain and Vi-Comp II Alloy.

Adverse Reactions

- Exposure to alloy dust or fumes may cause eye irritation and/or respiratory complications.
- Certain components of alloys are potential carcinogens. See Warnings.
- Ingestion of gallium may cause bone marrow depression.

STEP BY STEP INSTRUCTIONS

WAX-UP AND SPRUING:

Use a die spacer and wax copings to a minimum of 0.4mm. Indirect spruing for multiple units and direct for single units: use 10 gauge (2.6mm) sprue for copings; 8 gauge (3.2mm) sprues for bulky crowns and pontics. Sprues should be 6.5 to 9mm in length. Use 5 gauge (4.5mm) runner bars. The space between the runner bar and pattern should be 2-4 mm.

INVESTMENT:

Use a high heat phosphate investment such as Vi-Comp® Investment (ReOrder# N251500).

BURNOUT INSTRUCTIONS:

Bench set molds for one hour. Place molds in a cold furnace and heat molds to 1800°F (982°C) following the manufacturer’s recommendations. Hold for one hour.

TORCH CASTING:

Use a separate high heat crucible for melting Vi-Comp II Alloy. Use propane oxygen or oxygen acetylene multi-orifice torch. Wind casting machine one more turn than used for gold alloys. **Note: Do Not break oxide layer. The casting temperature is reached when the molten mass rolls or waves from pressure of flame.**

INDUCTION CASTING:

DO NOT USE CARBON LINER IN A CRUCIBLE. A preheated crucible will minimize the amount of metal left in the crucible after casting.

Machine Settings:

Casting Temperature Setting: . 2750°F (1510°C)
Soak Timer:..... 15 seconds
Pyrometer position: . Back (focused on the alloy)
Acceleration reducer:..... 10-30

SURFACE PREPARATION:

Grind all alloy surfaces, which are to receive porcelain, using aluminum oxide stones. Blast the surface using non-recycled 125 micron dental grade aluminum oxide, max 3 bar (44 psi). The framework should be steam cleaned or ultrasonically cleaned in distilled water. (DO NOT use soap, acids or solvents).

PREOXIDATION:

Fire under vacuum 1400°F (760°C) to 1800°F (982°C) at 90°F (50°C) per minute with no soak. After the units have cooled, remove the oxide by blasting using non-recycled 125 micron dental grade aluminum oxide, max 3 bar (44 psi). Clean the units as described above.

OPAQUE AND BODY PORCELAIN BUILD-UP:

Apply and fire opaque and body porcelains per manufacturer’s instructions.

SOLDERING (Recommended Solders & Fluxes):
Pre-Soldering:Use a gas-oxygen torch.

Pre-Solder:ReOrder# N764300
While Bake-On Pre-Solder

Pre-Soldering Flux: .ReOrder# N200122
Pre-Porcelain Soldering Flux
High Fusing

Post Solder:ReOrder# NO86100
.583 Fine Low Fusing
Gold Solder

Post Soldering Flux:ReOrder# N200224
Post-Porcelain Soldering Flux

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