**Indication:** Brazing material for dental alloys.  
**Contraindication:** If patients are known to be allergic to any of the components, the material should not be applied.  
**User tip for post-soldering (furnace soldering #):**  
 Correctly position the parts to be soldered on the master model. Fix the parts to be soldered with a synthetic that can be fully burnt out (alternatively, if this is not available, use a fully combustible adhesive wax). *see* ①  
 Make a durable soldering jig using a reliable soldering investment (or a reliable casting investment mixed with distilled water). NB: Ceramic parts must NOT come into contact with the investment. *see* ②  
 Remove the fixing material after the investment has hardened. *see* ③  
 Dry the (soldering) investment at a temperature no higher than 250°C to prevent unnecessary oxidation of the parts to be soldered. *see* ④  
 Apply soldering flux into the soldering joint and onto the solder.  
 - Use Elephant Paste Flux with **precious** metal combinations and Elephant Fluorflux for **non precious** metal combinations. *see* ⑤  
 Apply a few pieces of solder to the soldering joint and cover with soldering flux. Programme the porcelain furnace and place the object in the furnace using a firing tray. *see* ⑥ (See below for a programming example for soldering with a 750°C solder)

Start Temp	Hi T*	Dry	Hold	Heat	Vac.*	Rate of increase	Cool
500°C	795°C (=T+45°C)	2 - 4 min.	1.5 - 3 min.	1 min.	yes, (500°-795°C)	60°C / min.	0 min.

\* Adjust the temperature (T) for the other postsolders, T + 45°C.  
 Allow enough time for cooling the object and only then remove the investment soldering jig. *see* ⑦  
 Flux and investment residues are easily removed by sandblasting and pickling. Finish the solder area with the usual finishing tools and materials. *see* ⑧

**Please note:**

- Use the correct soldering flux.
- Do not melt the soldering wire into a spherical shape beforehand. Use pure soldering material.
- Avoid overheating the solder at all times.
- Allow enough time for cooling the preparation.
- In the case of furnace soldering with precious metals, it is advisable to place a block of active carbon near the soldering place to achieve an additional oxide reduction.

# For flame soldering: see the directions for use on the pre-solder packaging or go to [www.elephant-dental.com/downloads](http://www.elephant-dental.com/downloads).

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