Working procedure for CERA AURUM

- 1. Indication.** Light yellow dental alloy for metal-ceramic, inlays, crowns, bridges and telescopic technique.
- 2. Contraindication.** If patients are known to be allergic to any of the components, the material should not be applied.
- 3. Sprueing.** According to the indirect method: 3 mm \varnothing for crowns, 3.5 mm \varnothing for pontics, runner bar 4-4.5 mm \varnothing .
- 4. Investing.** Invest in a reliable phosphate-bonded investment, such as Carrara Universal DL Investment.
- 5. Burn out.** Preheat at 800°C/1470°F for 30-120 min. depending on the size of the muffle ring.
- 6. Casting.** Casting temperature: 1200°C/2192°F in a graphite crucible; or with the flame in a ceramic crucible with for example Elephant Proflux. Use at least 1/3 new metal at each casting. Clean the metal well before re-use, for example by means of sandblasting and steam-cleaning.
- 7. Finishing.** Use carbide burs and/or ceramic bonded stones. The final trimming must be carried out in one direction only by means of fine crossteethed carbide burs! Sandblast with 50-125 micron aluminumoxide (Alublast), (pressure max. 3 bar).
- 8. Cleaning.** With distilled water ultrasonic or by means of a steam cleaner.
- 9. Oxidation.** 3 minutes at 860°C/1580°F without vacuum.
- 10. Bake.** With Carrara ceramic, cooling stage; slow, 5-7 minutes.
- 11. Pre-soldering.** Carrara 950°C/1740°F with Elephant Paste Flux.
- 12. Post-soldering.** Carrara 750°C/1380°F with Elephant Paste Flux.

For a reliable quality we advise to use Cera Aurum in combination with Carrara ceramic.

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