Working procedure for ORION YELLOW

- 1. Indication.** Yellow dental alloy for metal-ceramic, inlays, crowns and bridges.
- 2. Contraindication.** If patients are known to be allergic to any of the components, the material should not be applied. Bridges which exceed 6 units.
- 3. Sprueing.** According to the indirect method: 3 mm \varnothing for crowns, 3,5 pontics, runner bar 4-4,5 mm \varnothing .
- 4. Investing.** Invest in a reliable phosphate-bonded investment, such as Carrara Universal DL Investment.
- 5. Burn out.** Preheat at 800°C/1470°F for 30-120 min. depending on the size of the muffle ring.
- 6. Casting.** Casting temperature is 1200°C/2190°F in graphite crucible; or with the flame in a ceramic crucible with Elephant Proflux. Use at least 1/3 new metal at each casting. Clean the metal well before re-use, for example by means of sandblasting and steam-cleaning.
- 7. Finishing.** Use carbide burs and/or ceramic bonded stones. The final trimming must be carried out in one direction only by means of fine cross-teethed carbide burs! Sandblast with 50-125 micron aluminium oxide (Alublast), (pressure max. 2 bar).
- 8. Cleaning.** With distilled water ultrasonic or by means of a steam cleaner.
- 9a. Oxidation.** 3 minutes at 900°C/1650°F, without vacuum.
- 9b. Oxidationreduction.** Sandblast with aluminum oxide 125 μm (pressure max. 2 bar) or pickle for max. 1 min. in clean pickling agent. Afterwards clean again (see point 8).
- 10. Bake.** With a standard fusing ceramic, such as Antagon Interaction. Baking temperature may not exceed 920°C/1688°F. If the porcelain manufacturer prescribes a baking temperature higher than 920°C/1688°F, but no higher than 930°C/1705°C, the baking time at 920°C/1688°F should be extended by one minute for every additional 10°C/50°F. With Antagon Interaction: 5-7 minutes cooling down.
- 11. Pre-soldering.** Carrara 950°C/1740°F with Elephant Paste Flux (no large surfaces).
- 12. Post-soldering.** Elephant II PdF 800°C/1470°F or III (PdF) 750°C/1380°F with Elephant Paste Flux.

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