## User instructions

# Vi-Comp NE bonder

A supplementary ceramic product for the Antagon Interaction ceramic system, is used to balance non-precious metal alloys and the ceramic veneering.

Made in Holland

Rx only  
ISO 9693 + ISO 6872

Verlengde Lageweg 10  
1628 PM Hoorn, The Netherlands  
Tel: +31 229 25 90 00  
Fax: +31 229 25 90 99  
E-mail: info@elephant.nl  
www.elephant-dental.com

## Elephant Dental B.V.

### Indication:

- Vi-Comp NE bonder, a supplementary ceramic product for the Antagon Interaction ceramic system, is used to balance non-precious metal alloys and the ceramic veneering.
- For ceramic alloys with a CTE of 13.8 – 14.9  $\mu\text{m}/\text{m.K}$  (25 - 500°C).
- For optimum results, we recommend the use of 100% new NP alloys.

### Contraindication:

- All other alloys not listed under 'Indication'.
- If patients are known to be allergic to one of the components, the material should not be applied.

### Preparation of the sub-structure:

- Finish the sub-structure according to criteria provided by the alloy manufacturer and avoid any sharp corners and edges.
- Clean the metal surface by blasting with 125  $\mu\text{m}$  aluminium oxide at a pressure of 3 - 4 bar.
- Unless expressly indicated otherwise in the instructions, an oxidation firing phase is not required and should only be used to check the sub-structure.

### Application and processing:

- Mix Vi-Comp NE bonder with Superwet or distilled water into a thin creamy substance.
- Apply the mixture to the metal in a non-masking layer (Washbrand)
- Fire the bonder at 980°C to an equally uniform layer according to the following table:

Start temp.	Dry	Vacuum start	Rate of heat increase	Final temp.	Hold time (with vacuum)	Hold time (without vacuum)	Cooling down
400°C	4 min.	400°C	60°C/min.	980°C	1 min.	1 min.	-

- After this firing phase, the bonder is uniformly sintered onto the sub-structure, which ensures equal vitrification of the metal surface. This creates the possibility of achieving reliable adhesion between an NE ceramic alloy and the ceramic facing.
- The values stated are target values and should be used exclusively as guidelines. Varying firing results are possible depending on the capacity, manufacturer and age of the ceramic oven being used.

### Ceramic processing:

- Follow this bonder firing with the Antagon Interaction instructions from step 2 of the process onwards (opaque firing procedure).
- To prevent unwanted ceramic discolouration, avoid direct contact between the ceramic facing and the alloy.

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