User instructions

Gebrauchsanleitung • Gebruiksaanwijzing • Instructions d'utilisation
• Istruzioni per l'uso • Instrucciones para el uso

ISO 22674

WIPTAM®

	REF.
round	7600141070
rund	7600141090
rond	7600141100
rond	7600141110
fondo	7600141120
redondo	7600141130
	7600141150

	REF.
half-round	0,75 x 1,50 7600142075
halbrund	0,87 x 1,75 7600142087
half-rond	1,00 x 2,00 7600142100
semi-rond	
media caña	

Co=40.0, Cr=20.0, Ni=16.0, Fe=14.0, Mo=7.25,

Mn=2.0, Rest Si

Dichte, S.G. 8.3 g/cm³

Vickershärte, HV10 375

Bruchdehnung, Elong. 20%

Rp-0.2 820 MPa

F max. 1.100 MPa

E-Modul 230 GPa

Ni

Rx only

Cobalt-based dental wire for clasps, retainers and friction-fit pins.

WIPTAM is a colour-fast and mouth-resistant wire whose specific characteristics provide considerable advantages when fabricating clasps, retainers and friction-fit pins. When WIPTAM wires are supplied medium hard (F max. 1,000 - 1,200 MPa). When WIPTAM wires are cold-worked, e.g. during bending they reach a spring hardness.

Unlike FeCrNi clasp wires, no difficulties will occur when soldering WIPTAM wire. The hardness that is lost due to the heat applied during soldering is regained during mechanical post-treatment (bending, trimming or polishing).

WIPTAM wire is also successfully used for orthodontic objects that require soldering.

Indication:

WIPTAM is a colour-fast and mouth-resistant cobalt-based wire whose specific characteristics provide considerable advantages when fabricating clasps, retainers and friction-fit pins.

Contraindications:

WIPTAM wire contains nickel. It should not be used by or for people with a known Nickel allergy, or with an allergy or sensitivity to one of the other components. Consult the MSDS for more information.

WIPTAM wire is not suitable for casting channels.

WIPTAM wire cannot be used for casting.

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Adverse effects:

In exceptional cases, certain components of the alloy may cause an allergic reaction or sick feeling due to electrochemical processes.

Interactions:

In exceptional cases, occlusal or approximal contact with different alloys may cause a sick feeling due to electrochemical processes.

Package:

The package contains eight straightened 25-cm wire rods.

Using WIPTAM:

To ensure the stability of the clasp basis, always use wire that is thicker than the clasp to be made. Before soldering, bend the clasp with smooth-tipped pliers. Sharp-edged tools or pliers with rough beaks may create grooves, which are potential weak spots. White gold solder or WITEX CoCr solder produce the best soldering results. After soldering, adapt the clasp optimally along the contours of the element. Cut the wire to the desired length and taper the end of the clasp. Bending, trimming and polishing WIPTAM wire changes its material structure, converting it from medium hard to springy hard. The greater the force exerted to bend and fit the wire, the greater the increase in wire hardness. If the clasp has become too hard at one particular spot due to intensive mechanical treatment, this specific spot can be heated to make it flexible again.

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Elephant Dental B.V.