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## SILFLEX ORANGE SPEED

Addition-polymerising duplication silicone for high-precision duplication moulds.

**Indication:** The two-component addition-polymerising duplication silicone SILFLEX ORANGE SPEED is the ideal material for producing accurate model reproductions for complicated high-precision dental work.

**Contra-indication:** If users are known to be allergic to any of the components, the material should not be used.

### Product data

Composition	Polymethylhydrogen-siloxanes
Colour	orange
Hardness	15 shore A
Mixing ratio	1 : 1
Working time (at 22°C room temperature)	approx. 3 min.
Model removable from mould at 22°C	after 8 min.
Shrinkage after 7 days	0.05%
Elongation at break (DIN 53 504, standard ring)	approx. 320%

### Processing

#### 1. Models

Prepare the (dry) plaster model in the usual way for duplication. Place the model in the duplicating flask and close. Make sure that the flask is well sealed and the models to be duplicated are at room temperature.

#### 2. Mixing

Using a spatula, mix the silicone components A (white) and B (orange) in a ratio of 1:1 (weight) for approx. 30 seconds until the mixture has a uniform colour, indicative of a homogenous mixture.

Tip: Use a silicone mixer to achieve a homogenous 1:1 mixture. This way, measuring errors are avoided and important material savings made.

#### 3. Casting

Fill the flask holding the dry plaster model with the mixed SILFLEX ORANGE SPEED duplicating silicone. To avoid air inclusions, ensure that the silicone compound is poured in a thin stream from a height of 30 to 40 cm onto a spot in the centre of the model so that it can flow evenly out from there. Make sure that the upper surface of the model is covered with a layer of silicone at least 5 mm thick. Do not bring the duplicating silicone into contact with hard rubber or plasticine as it cannot set at points of contact with these materials.

#### 4. Hardening

At room temperature, the silicone mould sets in 8 minutes. Any trapped micro-bubbles, which may form if the silicone compound is cast too quickly, will not be a problem if hardening takes place under pressure (3-4 bar, max. 30°C).

Note: All operations must be carried out under identical conditions. If the duplicating mould is produced under pressure, the duplicate model must harden under pressure as well.

#### 5. Removing the model from the silicone mould

The model can be taken from the silicone mould after a hardening period of 8 minutes (e.g. using compressed air).

#### 6. Fabricating the duplicate model

Spray the surface of the silicone mould with a surface tension reducer. Allow to dry properly before filling the mould with the model compound. All usual plasters, investments and modelling resins can be used to make models in SILFLEX ORANGE SPEED duplicating moulds. Keep the silicone mould in the flask until the model is fully hardened.

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