Indication: Brazing material for dental alloys.

Contraindication: If patients are known to be allergic to any of the components, the material should not be applied.

User tip for pre-soldering:

Correctly position the parts to be soldered on the master model. Fix the parts to be soldered with a synthetic that can be fully burnt out (alternatively, if this is not available, use a fully combustible adhesive wax). *see 1*

Carefully remove the preparation from the model without deforming it. Make a durable soldering joint using a reliable soldering investment (or a reliable casting investment mixed with distilled water). *see 2*

Remove the fixing material after the investment has hardened. *see 3*

Dry the (soldering) investment at a temperature no higher than 250°C to prevent unnecessary oxidation of the parts to be soldered. *see 4*

Apply soldering flux into the soldering joint and onto the solder.

- Use Elephant Paste Flux with **precious** metal combinations and Elephant Fluorflux for **non precious** metal combinations.

First heat the object evenly with the soldering flame. Once the object begins to glow evenly, the solder area itself can be heated further. *see 5*

Heat the solder area in such a way that the solder flows under the influence of the heat of the object. Because solder has the tendency to flow towards the heat source, hold the flame on the side opposite the side to which the solder is being applied. *see 6*

Once the soldering area has flowed sufficiently, remove the flame from the object. Allow enough time for cooling the object and only then remove the investment soldering jig. *see 7*

Flux and investment residues are easily removed by sandblasting and pickling. Finish the solder area with the usual finishing tools and materials. *see 8*

Please note:

- Use the correct soldering flux.
- Do not melt the soldering wire into a spherical shape beforehand. Use pure soldering material.
- Avoid overheating the solder at all times.
- Allow enough time for cooling the object.

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