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# R<sup>x</sup> only CE

## Warning: This product contains substance(s) known to the State of California to cause cancer, or other reproductive harm. (Calif. Prop. 65)

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Warning: This product contains substance(s) known to the State of California to cause cancer, or other reproductive harm. (Calif. Prop. 65)

Warning: This product contains substance(s) known to the State of California to cause cancer, or other reproductive harm. (Calif. Prop. 65)

	Manufacturer		EC REP
DENTSPLY International Inc.	Elephant Dental B.V.	DENTSPLY AUSTENAL	Verlengde Logweg 10
York, PA 17405-0872	1628 PM Hoorn	1-800-621-0381	The Netherlands
www.austenal.dentsply.com	+31 229 25 90 00	www.dentsply.com	+229 259 0000
			www.elephant-dental.com

Vitallium<sup>®</sup> Trademark is licensed to DENTSPLY Austenal by SMD Corporation. Perflex<sup>®</sup> and Silflex<sup>®</sup> are registered trademarks of DENTSPLY Austenal.

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
**Indications For Use**

Vitallium<sup>2</sup> Alloy is a cobalt chromium molybdenum based casting alloy designed for use when casting removable partial denture frameworks.

**Contraindications**

- Vitallium<sup>2</sup> Alloy is contraindicated for patients and users with a history of sensitivity to cobalt, chromium, and molybdenum.

**Warnings**

- Prolonged exposure to alloy dust and/or fumes may lead to lung irritation, and/or pulmonary complications. Use appropriate engineering controls to limit exposure. For excessive inhalation of dust or fumes, seek medical advice.
- This alloy contains chromium. Some compounds of these elements are potential carcinogens. Advise your health care provider of exposure to these elements.
- This alloy contains elements that are known to the State of California to be carcinogenic.

**Precautions**

- When melting, grinding, or polishing alloys, use adequate ventilation, vacuum systems, protective eyewear, and protective masks and clothing.

**Adverse Reactions**

- Exposure to alloy dust or fumes may cause eye irritation and/or respiratory complications.
- Certain components of alloys are potential carcinogens. See Warnings.

## STEP BY STEP INSTRUCTIONS

**PROCESSING INSTRUCTIONS**

**Investment Recommendations:**

Use a DENTSPLY Austenal’s Ethyl Silicate Bonded Investment System.

ReOrder#

VR™ Investment ..... NO22065

VR Binder ..... NO23001

VR Binder Thinner ..... NO24001

**Duplicating Material Recommendations:**

Use **Perflex®** Reversible Hydrocolloid Duplicating Material (ReOrder# NO21055) with **VR Investment**.

**Silflex®** III Polyvinyl Siloxane Duplicating Material (ReOrder# NO21501) is recommended for use for model duplication for cases with precision attachments and for combination cases.

**NOTE: Follow the instructions for the investment and duplicating material that you are using for refractory model and investing steps.**

### EQUIPMENT

**BURNOUT FURNACE REQUIREMENTS**

**For VR Investment**

High temperature gas fired burnout furnace capable of attaining 2150°F (1177°C).

**CASTING MACHINES**

DENTSPLY Austenal ECM IV Casting Machine or a Torch Casting Machine - i.e. Oxygen acetylene multi-orifice torch & spring activated centrifugal casting arm.

**BURNOUT INSTRUCTIONS**

Preheat the investment molds and a clean casting crucible following the recommended burnout cycle for the investment being used. *Make sure that crucibles are only used for this alloy and that used crucibles are clean of slag and any remaining metal.*

**CASTING INSTRUCTIONS**

**Virgin Metal/Revert Ratio:** To maintain the physical and mechanical properties of this alloy, it is recommended that at least 50% new metal be used for each casting. Revert (buttons) must be sandblasted and cleaned prior to reuse.

By adhering to this recommendation, the chemical composition of the alloy will be maintained and excess oxide or slag build up in the alloy (a major cause of pits in castings) will be avoided.

**DENTSPLY AUSTENAL ECM IV CASTING MACHINE SETTINGS**

Casting Temperature Setting: ..... 2820F Vitallium Scale
Soak Timer: ..... 5 seconds
Melt Selector Switch: ..... On High

**ACCELERATION REDUCER SETTING:**

Small Castings - 50 Large Casting - 75 to 100 RPM: 380-400

*Adjust centrifugal arm counter weight to the proper position to provide for a smooth and balanced casting operation.*

Follow the instructions for the machine that you are using. Melt and cast the alloy into the preheated investment mold and allow the molds to bench cool before further processing.

DO NOT FORCE COOL OR QUENCH THE MOLDS AS IT MAY CAUSE WARPAGE OF THE CASTINGS & UNDESIRABLE CHANGES IN THE MECHANICAL PROPERTIES OF THE ALLOY.

**TORCH CASTING INSTRUCTIONS**

**Warning:** *Use extreme caution during torch melting operations. Follow the instructions and warnings of the torch manufacturer and wear eye protection designed for use when operating an oxygen acetylene torch.*

Use 3-5 turns on the casting arm depending on the tension of the spring. Set the counter weights at the correct position according to the size of the case to be cast.

Use a multi-orifice oxygen acetylene torch pressure regulators and flow meters.

Place the preheated mold and crucible in the casting machine. Place the required amount of alloy in the crucible.

Light the torch and with a circular motion move the flame slowly over the metal. Continue this motion to uniformly heat and melt the alloy. When the alloy begins to melt, the ingots and buttons will collapse. Add flux and continue to uniformly heat the alloy until the molten alloy can be seen to slump and move under the pressure of the flame. DO NOT overheat the alloy.

**Note: The oxide film should not break.**

Release the casting arm and allow the arm to coast to a complete stop. Remove the mold and allow it to bench cool before further processing.

DO NOT FORCE COOL OR QUENCH THE MOLDS AS IT MAY CAUSE WARPAGE OF THE CASTINGS & UNDESIRABLE CHANGES IN THE MECHANICAL PROPERTIES OF THE ALLOY.

**INSTRUCTIONS FOR CLEANING PARTIAL DENTURES**

Vitallium<sup>2</sup> Alloy dentures can be cleaned using commercial denture cleansers or denture pastes which indicate on their labels that they are suitable for cleaning “bridges, orthodontic appliances and/or partial dentures.” The use of cleansers having high hypochlorite contents and the use of household bleach is not recommended as they may cause tarnish or pitting of the framework.

**INSTRUCTIONS FOR THE REPAIR OF VITALLIUM<sup>2</sup> ALLOY FRAMEWORK**

Frameworks can be repaired using Vitallium Welding Rod (ReOrder# NO41200 & 01) or Nickel Free Welding Rod (ReOrder# NO41005) or White Gold Solder (ReOrder# NO85100 or NO85200).

**PHYSICAL & MECHANICAL PROPERTIES**
0.2% Yield Strength ..... 80,000 psi (551 MPa)
Ultimate Tensile Strength.....109,000 psi (751 MPa)
Modulus of Elasticity .....29,000 kpsi (200 GPa)
Percent Elongation ..... 9.6%
Vickers Hardness ..... 395 HV5
Density ..... 8.3 g/cc
Melting Range .....2375-2500°F (1300-1370 °C)

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