User instructions

CHROMODUR®

Cobalt based dental casting alloy for removable restorations.

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Indications:
Type V, dental casting alloy of cobalt, chromium and molybdenum for model casting whose excellent, finely-balanced properties considerably exceed the values specified in ISO 22674 and international standards. The use of materials largely free from impurities and special melting and manufacturing processes impart extreme purity to CHROMODUR. It contains neither nickel nor beryllium as alloying elements and displays excellent tissue compatibility. CHROMODUR alloys are very easy to work with.

Contraindications:
If components of the alloy are not tolerated or cause an allergic reaction, use must be terminated.

Adverse effects:
In exceptional cases, certain components of the alloy may cause an allergic reaction or sick feeling due to electrochemical processes.

Interactions:
In exceptional cases, occlusal or approximal contact with different alloys may cause a sick feeling due to electrochemical processes.

Processing instructions

1. Sprues
Apply 2 to 4 sprues to each wax pattern, depending on the size of the object.
Sprue diameter: 3-3.5 mm.
Connect at the pattern's thickest point expanding in delta fashion.
Arrange sprues in line with the direction of metal flow. Keep sprues short to minimize distance between cone and cast. Attach sprue former at least 1 cm above pattern.

2. Painting with liquid investment
Apply a thin, even coat of liquid investment to the pattern, avoiding bubbles.
Invest model immediately, do not allow liquid investment to dry.

3. Investing
PH2 is ideally suited to reach optimal fit. Pour into casting ring, avoiding bubbles.
Mixture ratio powder : liquid
Models [with duplicating compound] 100 g : 13 ml
[with duplicating silicone] 100 g : 16 ml
Moulds 100 g : 16 ml

hardening time: 30-45 min.
(see Processing instructions PH2)
Alternatively to PH2, we also offers the PH3 binder investment material for the CHROMODUR system.

4. Preheating
Preheat casting rings in preheating furnace according to PH2 processing instructions:
300°C (572 °F) holding time 30 min.
600°C (1,112 °F) holding time 30 min.
1,000°C (1,832 °F) holding time 30-45 min.
For plates and delicate designs preheat at 1,050°C (1,922 °F) holding time 30-45 min.

5. Required metal quantities
Simple bar: 3 ingots = 25 g
Skeletal plate: 3-4 ingots = 25-34 g
Complete plate: 5-6 ingots = 42-51 g

6. Casting
a) Motorized or spring-loaded centrifugal caster for open-flame melting
Ignite the acetylene-oxygen or propane-oxygen gas mixture. Fully open both valves and adjust the flame by throttling the intake of acetylene or propane such that blue cones approx. 3 mm long are visible at the holes of the multi-orifice head. Pressure values for the flame: acetylene = 0.7 bar, oxygen = 1.9 bar or propane = 0.3 bar, oxygen = 1.5 bar. This gives a neutral flame which is absolutely essential for melting CHROMODUR. Place CHROMODUR ingots into preheated crucible. Heat ingots evenly until they form a molten bath by keeping the torch about 40 mm above the crucible whilst making rotating movements. As soon as the oxide shade breaks, initiate centrifugal movement.

b) Atmospheric high-frequency melting and centrifugal casting unit
Do not use graphite crucibles.
As soon as ingots are molten and shaded areas have disappeared from the melt, wait 4 sec and then initiate centrifugal casting.

c) Vacuum high-frequency melting and centrifugal casting unit
This unit can be used for melting and casting all CHROMODUR alloys under vacuum conditions. As soon as ingots are molten and shaded areas have disappeared from the melt, wait 4 sec and then initiate centrifugal casting.

See www.elephant-dental.com for material safety data sheets (MSDS).

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