DENTSPLY
PROSTHETICS

EVER *soft*®

SOFT DENTURE LINER

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Indications

EverSoft® is a two component (powder and liquid) plasticized methacrylate indicated for use in stabilizing and relining complete and partial dentures, cushioning sensitive tissues and retaining overdentures. The compliance and resiliency of EverSoft can be varied by altering the powder/liquid ratio. EverSoft mixed using a powder/liquid (vol/vol) ratio of 2:1 results in a very compliant lining layer suitable for post-surgical prosthetics or as a cushion against sensitive tissues. EverSoft mixed using a powder/liquid ratio of 2 1/2:1 results in a firmer lining which effectively stabilizes partial and complete dentures while still providing the advantages of a soft-liner for long-term use. EverSoft mixed using a powder-liquid ratio of 3:1 results in a lining capable of retaining overdentures onto natural teeth, ball abutments or bar connectors. EverSoft forms a durable bond with conventional methacrylate resin denture base formulations.

Contraindications:

None known.

Warnings:

EverSoft contains methacrylate resins which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical attention. EverSoft liquid contains ethyl alcohol and should be kept away from sources of ignition. Avoid contact with eyes. In case of accidental contact, wash thoroughly with water and seek medical attention.

Keep containers tightly closed when not in use.

Precautions:

Work in a well ventilated area when applying EverSoft Sealer and DO NOT place prosthesis in patient's mouth until sealer is completely dry and prosthesis has been rinsed with water. Consult the Material Safety Data Sheet (MSDS) for more information.

It's longevity of service is reduced when used as a long-term soft liner against sharp, knife edge ridges. It can be used for application over knife like ridges, but use a powder/liquid ratio of 3:1 for more resistance to hairline cracking. Adhesion of EverSoft to each and every denture base resin formulation has not been tested. Carefully evaluate adhesion of EverSoft to denture base formulations which do not contain polymethylmethacrylate.

Adverse Reactions:

Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.

LABORATORY TECHNIQUE

Density	Ratio (P:L)	Powder	Liquid	Applications
Softest	2:1	8cc	4cc	<ul style="list-style-type: none"> • Post-surgical prosthetics • Cushion for sensitive tissue
Standard Lining	2½ :1	10cc	4cc	<ul style="list-style-type: none"> • Complete denture relines • Partial denture stabilization • Overdenture retained on healing caps
Firm	3:1	12cc	4cc	<ul style="list-style-type: none"> • Overdenture retained by ball abutments, bar connectors or teeth.

Preparation:

1. Upon receipt of reline impression, pour a stone model. Flask or jig the model with the impression in place following standard reline procedures.
2. After the flask has set for one hour, place it in a boilout tank for approximately 15 minutes. Remove it and separate the upper and lower portions of the flask.
3. Gently scrub the lower half of the flask (model half) with a mild detergent and rinse thoroughly with clean, hot water. Post dam, if necessary; apply separator; allow the separator to dry.
4. Prepare the denture base for the soft liner by relieving the acrylic tissue-bearing surface by 2mm or more. Liner can extend above the peripheral roll to any length.
5. Eliminate severe undercuts and roughen areas of denture base to be relined.
6. Gently scrub the upper half of the flask (denture half) and rinse with hot water. Apply separator to stone or plaster; let dry.

Measuring & Mixing:

7. To measure: use the graduated glass vials marked "Powder" and "Liquid." You can control the liner's softness

by modifying the powder/liquid ratio as shown in the chart above.

Note: At a 2:1 ratio, liner is the softest. A light powdering of liner's surface will absorb any excess liquid.

8. Pour powder (all at once) into the liquid in mixing cup; *Stir very slowly* (to avoid the introduction of bubbles) until all of the powder particles are totally moistened. The mixture consistency instantly turns "sticky" and "honey-like."

Application:

9. *Immediately*, spread the mixture onto the totally dry denture base surface.
10. Close the flask or jig for curing. In jig, use fingertip coated with the mixing liquid to taper and smooth any excess material 3mm onto labial and buccal surfaces.

Curing:

11. Place flask onto a press; slowly apply pressure until the flask is completely closed.

Note: Do not trial pack the soft liner.

- 12a. In a flask: Place flask in a curing unit filled with boiling water (212°F/100°C) for 15 minutes.

- 12b. In a jig: Place jig in a curing unit or in a pressure pot at (20psi/1.4bar) in steaming hot water (between 122°F/50°C to 165°F/74°C) for 10 minutes.

Finishing:

- Remove cured denture from flask or jig; chill it a few minutes in an ice water bath to temporarily harden the soft liner for easy trimming.
- Remove excess flash with sharp scissors, scalpel or suitable rotary instrument (i.e. acrylic bur, rubber point or finishing stone). Feather taper liner 3mm onto labial and buccal surfaces.
- If desired, finish with a wet rag wheel and pumice to further smooth the surface. Do not high-shine the liner using a buff wheel as EverSoft Sealer provides a high-gloss finish.

Sealing:

- Dry the liner surface using air or warm air blower to remove all surface moisture. Apply a generous coat of EverSoft Sealer over the totally dry soft liner and air dry for 2+ minutes. Repeat procedure adding a second coating. A third coat of Sealer can be applied on non-retentive fitting surfaces.

Additional Information:

Note: Provide Patient Care Card to instruct care & cleaning of the soft liner.

Note: To adjust, refit or repair: roughen the surface, apply a fresh mix of soft liner, follow directions above for Application, Curing, Finishing and Sealing.

Note: The Powder is available in a "clear" shade which can

be used to retain transparent acrylic mouth guards, bite plates and splints, to protect veneers, and as a gasket connector around teeth.

Lining a New Denture:

Preparation:

- Invest denture, wax-up and boil-out in the usual manner. Apply separator to denture half of mould. Allow the separator to dry.
- Form a spacer on the model from 1mm to 2+mm in thickness. Allow spacer to overlap the peripheral roll 3mm onto the labial and buccal surfaces. Soft liner can extend above the peripheral roll to any length.
- With spacer in place, trial pack the acrylic with dry sheets. After final packing of acrylic, allow 20 minutes to "snap set" the acrylic while flask is in the press.
- Open flask, remove spacer and clean model. Apply separator onto the model-side stone surfaces. Allow the separator to dry.

Measuring & Mixing Application:

(Follow instructions on page 2.)

Curing:

Cure per manufacturer's instruction for the "new" acrylic resin being used.

Finishing and Sealing:

(Follow instructions on page 3.)

CHAIRSIDE TECHNIQUE

Density	Ratio (P:L)	Powder	Liquid	Applications
Softest	2:1	8cc	4cc	<ul style="list-style-type: none"> • Post-surgical prosthetics • Cushion for sensitive tissue
Standard Lining	2½ :1	10cc	4cc	<ul style="list-style-type: none"> • Complete denture relines • Partial denture stabilization • Overdenture retained on healing caps
Firm	3:1	12cc	4cc	<ul style="list-style-type: none"> • Overdenture retained by ball abutments, bar connectors or teeth.

Preparation:

1. Relieve all undercuts and remove 2mm or more of acrylic on the tissue side of the denture to make room for the soft liner.
2. Roughen the peripheral roll plus 3mm of the buccal and labial surfaces.
3. Clean and dry the areas where the soft liner is to be applied. Use petroleum jelly or a water soluble substitute to coat the interproximal surfaces that are not meant to be lined. This will ease finishing.

Measuring & Mixing:

4. To measure: Use the graduated glass vials marked "Powder" and "Liquid." You can control liner's softness by modifying the powder/liquid ratio as shown in the chart above.

Note: At a 2:1 ratio, liner is the softest. A light powdering of liner's surface will absorb any excess liquid.

5. Pour powder (all at once) into the liquid in mixing cup; Stir very slowly (to avoid the introduction of bubbles) until all of the powder particles are totally moistened. The mixture consistency instantly turns "sticky" and "honey-like."

Application:

6. *Immediately*, spatulate and spread the mixture onto the totally dry denture.
7. Apply lubricant to patient's lips or ask patient to moisten lips as the soft liner does not adhere to wet areas.

Note: For patients with sore spots or tissue/taste sensitivity, rinse the mouth for 20 seconds with an anesthetic mouth rinse.

8. Place denture in the mouth and instruct the patient to close gently into occlusion for one minute while performing muscle moulding at the same time, if desired.
9. Remove denture from the mouth. The soft liner should completely overlap the peripheral roll.
10. Use a fingertip coated with the mixing liquid to taper and smooth the excess material 3mm onto the labial and buccal surfaces; then cure immediately.

Curing:

11. Place the denture, reline side up, into a cup filled with steaming Hot Water (between 122°F/50°C to 165°F/74°C) and cure for 10 min. If available, place in a pressure pot;

cure in Hot Water (between 122°F/50°C to 165°F/74°C) at (20 psi/1.4 bar) for 10 min.

12. After curing, chill the denture in Cold Water or Ice Water or spray the liner's surface with a skin refrigerant to temporarily harden the soft liner for easy trimming.

Finishing & Sealing:

13. Remove excess flash with sharp scissors or scalpel. Feather taper the material onto the labial and buccal surfaces with an acrylic bur and smooth with a large rubber point or finishing stone.
14. If desired, finish with a wet rag wheel and pumice to further smooth the liner surface. Do not high-shine the soft liner.
15. Dry the liner surface using air or warm air blower to remove all surface moisture. Apply a generous coat of EverSoft Sealer over the totally dry soft liner and air dry for 2+ min. Repeat procedure adding a second coating. A third coat of EverSoft Sealer can be applied on non-retentive fitting surfaces.

DO NOT place denture in patient's mouth until EverSoft Sealer is completely dry and rinsed.

Patient Delivery:

Give patient a Patient Care Card to instruct him/her on care and cleaning of the soft liner.

- 1) Soak only 10 min. in denture cleaner.
- 2) Warn patient not to scrub the soft liner with a hard-bristle brush; and
- 3) Request annual re-sealing of soft liner and dental examination to assure proper fit.

Note: To adjust, refit or repair; roughen the surface, apply a fresh mix of soft liner, insert, cure, finish and re-seal.

Note: Liner will cure in 5 minutes in-the-mouth for thin, small area relines.

Note: The Powder is available in a "Clear" shade which can be used to retain transparent acrylic mouth guards, bite plates and splints, as a gasket connector around teeth, and to fill voids to eliminate food traps.

Implant Overdenture Technique:

This patented implant overdenture technique was developed to aid in the progressive healing process following the surgical placement of implants. It is well suited for this chairside application due to one's ability to vary the liner's degree of softness, (see mixing ratios and applications in chart on page 4.)

Post-Surgical Placement:

After a one to two week interval following the implant placement, line the tissue side of the overdenture for a soft, light load.

Note: Initially, use softest liner mix ratio of 2 parts powder to 1 part liquid.

After suture removal, block out suture line (indentation) with dry foil or suitable material and reline following the Chairside Technique. Make sure soft material does not interfere with healing site or tissue flap. This soft liner cushions tender areas and its sealer provides a non-absorbent, hygienic surface.

Healing Cap Exposure:

After an appropriate implant integration

phase, perform the “second stage” reline. This requires a new application of soft liner to the relieved tissue side to anchor over healing caps.

Note: Use the standard liner mix ratio of 2 1/2 parts powder to 1 part liquid.

This Healing Cap overdenture enables a “gradual healing load” with light functional stimulus to the implant. This thin area reline will cure in-the-mouth in 5+ minutes.

Long Term Retention:

After attachment of ball abutments or bar connectors, reline overdenture to form gasket over balls or bars “without” O-rings or clips.

Note: Use the firmest liner mix ratio of 3 parts powder to 1 part liquid.

The liner absorbs torque pressure and provides retention for long term service with easy chairside adjustment.

TECHNIQUE FOR CLEAR POWDER

Density	Ratio (P:L)	Powder	Liquid	Applications
Standard Lining	2 1/2 :1	10cc	4cc	Retain and cushion of clear acrylic appliances, such as night guards, bite plates, splints, TMJ & sleep apnea devices
Firm	3:1	12cc	4cc	<ul style="list-style-type: none"> • Gasket in a removable denture for retention on existing teeth • Fills voids beneath clasps and under retainers to eliminate food traps on acrylic prosthesis

Introduction

Long-lasting and easy to use, EverSoft Clear Powder possesses the same working characteristics as EverSoft Pink Powder. Once cured, the transparency of EverSoft Clear Powder makes it ideal for use with any clear acrylic resin appliance (e.g. night guards, bite plates) or as a gasket connector around teeth in a removable denture.

Soft Lining of Clear Acrylic Appliance

1. Process acrylic appliance in usual manner.
2. Relieve 2mm of acrylic from the occlusal surface of the appliance to make room for the soft lining.

Note: Maintain an acrylic border no thinner than 1mm to support soft lining.

3. Clean and thoroughly dry areas where soft liner is to be applied.
4. To place soft lining, refer to directions for “Laboratory Technique” or “Chairside Technique”.

Creating a Gasket in Removable Denture

1. Create a hole in denture with 2mm to 3mm of space around existing tooth or teeth to allow room for EverSoft gasket.
2. Roughen 3mm area surrounding the newly-created hole(s) on all surfaces where the soft liner will adhere.

3. To measure: Use the graduated glass vials marked "Powder" and "Liquid". For this application we recommend a 3:1 Powder-to-Liquid ratio.
4. Apply lubricant to all areas where you do not want the soft liner to adhere.
5. Follow Mixing & Measuring instructions in Directions for Use.
6. Spatulate mixture onto the tissue side of the totally dry denture beneath created hole(s).
7. Seat the denture over the existing tooth or teeth. The soft liner will flow out through the hole.
8. Use a fingertip coated with the mixing liquid to taper and smooth excess soft liner 3mm onto the labial and buccal surfaces.
9. Follow Curing, Finishing & Sealing instructions in Directions for Use.

Soft Lining of New Appliance

Follow Directions for Use, "Laboratory Technique," Lining New Denture, in order to place EverSoft Clear inside transparent acrylic appliances such as a night guard or bite plate. Make sure to leave 1+mm of acrylic border to support the soft material inside the arch.

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DENTSPLY International Inc.
Elephant Dental B.V.